

Leisure COLLECTION



shown in
Rosewood

Autumn



Cornflower



Spring



Rosewood



Chestnut



Burgundy



The comfort of the perfect sized chair has expanded into the Leisure collection! The UC332 with a well appointed, seams back and stylish armrest design is available in medium and large sizes for optimized comfort at an affordable price.

UC 332-M

UltraComfort[™]
AMERICA

Positioning you for life.

PROUDLY

MADE IN AMERICA
with Domestic & Imported Materials

Leisure COLLECTION

UC 332M

Model Number	UC332M
Size	Medium
Weight Capacity	375 LB
Recline Positions	3
Standard Fabrics	6
Chaise Seat	No
Back Type	Seam
Head Rest Cover	Yes
Arm Covers	Yes
Battery Back Up	Yes
Overall Width	32"
Overall Height	42"
Width Between Arms	21"
Floor to Top of Seat	18.5"
Seat Depth	21.5"
Seat to Top of Back	28"
Seat to Top of Arm	7"
Seat Angle Raised	16°
Floor to Top of Foot Rest Reclined	23"
Floor to Top of Back Reclined	33"
Back Angle Closed	110°
Back Angle Reclined	140°
Foot Rest Angle Closed	90°
Foot Rest Angle Opened	176°
Distance Required from Wall Reclined	17"
Extension Reclined from Back of Seat to End of Foot Rest	43"
Life Time Warranty, Chair Frame, Lift Frame, Recline Mechanism	Yes
Warranty on Electrical parts	2 YR
Pro Rated Years 4 Thru 7	No
Warranty on Mechanical Labor	1 YR
Warranty on Electrical Labor	1 YR
Warranty on Fabric and Foam	1 YR
Removable Back Option	Yes
Customizable *	No
Left Hand Control Option Available	Yes
Extra Pocket Option Available	Yes
4-Motor Heat & Massage Available	Yes
2-Motor Heat & Massage Available	No
Foot Rest Extension Option Available	Yes
Smart Tek System	Yes

* See Customize Chart for dimension limits
 Sizes may vary depending on fabric, filling material, upholstery and carpet thickness. All measurements taken on concrete floor using levels and metal rulers.



User-friendly feature of every UltraComfort lift and recline chair:

Easy to use battery backup system. Each UltraComfort lift and recliner has a battery backup system in the electrical transformer which operates on two 9 volt batteries. In the event of a power failure, the battery system will raise the chair from the reclined to the seated position. The batteries should be replaced after one use.



UltraComfort™
 AMERICA

Positioning you for life.

UltraComfort
 Phone: 866-793-5100
 Fax: 866-255-9962
 web: www.ultracomfort.net



Literature is current at time of printing. UltraComfort reserves the right to make changes to the product or literature at any time.



Complies with UL 73 and CSA C22.2 No. 68

Listing # E112417

FNAE 80-214A



SHT-UC332M.052013